

Anti-FILIP1L Antibody (3N116)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3N116
Purification:	Protein G purified

Applications

Western Blot

-Positive WB detected in: 293 whole cell lysate, MCF-7 whole cell lysate, PC-3 whole cell lysate, Hela whole cell lysate, HepG2 whole cell lysate

Verified Activity:	-All lanes: FILIP1L antibody at 1:1000
	-Secondary: Goat polyclonal to mouse IgG at 1/50000 dilution
	-Predicted band size: 102 kDa
	-Observed band size: 29, 80, 102, 130 KDa

-Exposure time: 5min

Application:	ELISA,WB
--------------	----------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Filamin A-interacting Protein 1-like Protein (1-200AA)
Antigen Species:	Human
Gene ID:	11259
Uniprot ID:	Q4L180
Biology Area:	Cardiovascular

Research Background

Acts as a regulator of the antiangiogenic activity on endothelial cells. When overexpressed in endothelial cells, leads to inhibition of cell proliferation and migration and an increase in apoptosis. Inhibits melanoma growth When expressed in tumor-associated vasculature.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481