

Anti-FLI1 Antibody (6L679)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6L679
Purification:	Affinity-chromatography

Applications

Western Blot	-Positive WB detected in: A549 whole cell lysate, U-251 whole cell lysate, Hela whole cell lysate
Verified Activity:	-All lanes: FLI1 antibody at 1:1000 -Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution -Predicted band size: 51, 44, 48, 30 kDa -Observed band size: 51 kDa
Application:	ELISA,WB
Recommended	WB:1:500-1:5000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human FLI1
Antigen Species:	Human
Gene ID:	2313
Uniprot ID:	Q01543
Synonyms:	Friend leukemia integration 1 transcription factor;Transcription factor ERGB;FLI 1;Proto-oncogene Fli-1
Biology Area:	Epigenetics and Nuclear Signaling, Cancer, Cell biology, Metabolism

Research Background

Sequence-specific transcriptional activator. Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481