

Anti-Ferritin light chain Antibody (6E431)

Product Details

Ig Type:	IgG2b
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6E431
Purification:	Protein G purified

Applications

Verified Activity:	<ol style="list-style-type: none">1. Immunohistochemical of paraffin-embedded human liver tissue using TMAH-00461 at dilution of 1:2002. Immunohistochemical of paraffin-embedded human brain tissue using TMAH-00461 at dilution of 1:2003. Immunohistochemical of paraffin-embedded human lymph node tissue using TMAH-00461 at dilution of 1:200
Application:	ELISA,IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Ferritin light chain
Antigen Species:	Human
Uniprot ID:	P02792
Synonyms:	NBIA3;MGC71996;LFTD;ferritin, light polypeptide
Biology Area:	Neuroscience

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481