

Anti-G4S linker Antibody-FITC (7T921)

Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	FITC
Clone:	7T921
Purification:	Protein G purified

Applications

Application:	FCM
Recommended	FCM:1:50-1:200.

Properties

Purity:	>95%
Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	G4S-BSA
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481