

## Anti-Growth Hormone Antibody (3P743)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3P743
Purification:	Protein G purified

### Applications

Application:	ELISA
--------------	-------

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human growth hormone
Antigen Species:	Human
Gene ID:	2688
Uniprot ID:	P01241
Synonyms:	IGHD1B;Growth hormone 1;GH1;hGH-N;GHN;GH-N;GH;GHB5
Biology Area:	Signal Transduction

### Research Background

Plays an important role in growth control. Its major role in stimulating body growth is to stimulate the liver and other tissues to secrete IGF-1. It stimulates both the differentiation and proliferation of myoblasts. It also stimulates amino acid uptake and protein synthesis in muscle and other tissues.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481