

Anti-GPRC5D Antibody (4T432)

Product Details

Ig Type:	Human IgG4(S228P)
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4T432
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Human GPRC5D with Anti-GPRC5D recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human GPRC5D at 5 µg/mL can bind Anti-GPRC5D recombinant antibody, the EC50 is 8.847-10.38 ng/mL.
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human GPRC5D Protein
Antigen Species:	Human
Gene ID:	55507
Uniprot ID:	Q9NZD1
Synonyms:	GPRC5D
Biology Area:	Immunology

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481