

## Anti-GPX1 Antibody (8R91)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8R91
Purification:	Affinity-chromatography

### Applications

Application:	ELISA
--------------	-------

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthetic peptide: Human GPX1
Antigen Species:	Human
Gene ID:	2876
Uniprot ID:	P07203
Biology Area:	Signal Transduction

### Research Background

Protects the hemoglobin in erythrocytes from oxidative breakdown. In platelets, plays a crucial role of glutathione peroxidase in the arachidonic acid metabolism.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481