

Anti-GRIA2/GRIA3 Antibody (1J893)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1J893
Purification:	Affinity-chromatography

Applications

Verified Activity:	Immunofluorescence staining of MCF-7 cells with TMAH-00514 at 1:56, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L).
Application:	ELISA,IF
Recommended	IF:1:20-1:200.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human GRIA2/GRIA3
Antigen Species:	Human
Uniprot ID:	P42262/P42263
Synonyms:	GRIA3;GLURC;AMPA-selective glutamate receptor 3;GluR-3;GluR-C;GluR-K3;Glutamate receptor ionotropic, AMPA 3;GluA3;Glutamate receptor 3;GLUR3
Biology Area:	Neuroscience

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481