

Anti-Hemoglobin Antibody (9J628)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9J628
Purification:	Protein G purified

Applications

Verified Activity:	1. Immunohistochemical of paraffin-embedded human spleen tissue using TMAH-00539 at dilution of 1:200
	2. Immunohistochemical of paraffin-embedded human placenta tissue using TMAH-00539 at dilution of 1:200
Application:	ELISA,IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Natural Human Hemoglobin
Antigen Species:	Human
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481