

Anti-Histidine-rich protein PFHRP-II Antibody (81865)

Product Details

Ig Type:	IgG1
Reactivity:	Plasmodium falciparum
Conjugation:	Unconjugated
Clone:	81865
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Histidine-rich Protein PFHRP-II for Plasmodium falciparum
Antigen Species:	Plasmodium falciparum
Uniprot ID:	P05227
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481