

## Anti-Histone H1.0 Antibody (4W4)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4W4
Purification:	Affinity-chromatography

### Applications

Application:	ELISA
--------------	-------

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthetic peptide
Antigen Species:	Human
Gene ID:	3005
Uniprot ID:	P07305
Synonyms:	H1FV;H1 histone family, member 0;SPG13;GROEL;HSP65;HSP60;HLD4;CPN60;H10;HSPD1;HuCHA60
Biology Area:	Epigenetics and Nuclear Signaling

### Research Background

Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The histones H1.0 are found in cells that are in terminal stages of differentiation or that have low rates of cell division.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481