

## Anti-Interferon alpha 2/IFNA2 Antibody (4U639)

### Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4U639
Purification:	Affinity-chromatography

### Applications

Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Mouse IFNA2 Protein
Antigen Species:	Mouse
Gene ID:	#N/A
Uniprot ID:	P01573
Synonyms:	LeIF A;IFN $\alpha$ 2;Interferon $\alpha$ -2;Interferon Alpha-2;IFN- $\alpha$ -2;IFNA2;Interferon Alpha-A;Interferon $\alpha$ -A;IFN-Alpha-2
Biology Area:	Immunology

### Research Background

#N/A

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481