

Anti-IFN gamma Antibody (4H542)

Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4H542
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human cervical using TMAH-00596 at dilution of 1:200
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Interferon gamma Protein
Antigen Species:	Human
Uniprot ID:	P01579
Synonyms:	IFN gamma;IFI;Interferon Gamma;IFG;IFN γ ;Interferon γ
Biology Area:	Immunology

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481