

Anti-IL-12RB1 Antibody (4C638)

Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4C638
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Human IL12RB1 with Anti-IL12RB1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized human IL12RB1 at 2 µg/mL can bind Anti-IL12RB1 recombinant antibody. The EC50 is 1.578-2.143 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human IL12RB1 Protein
Antigen Species:	Human
Uniprot ID:	P42701
Synonyms:	IL12RB;IL-12R-β1;interleukin 12 receptor, β1;interleukin 12 receptor, beta 1;IL-12R-BETA1; CD212;IMD30
Biology Area:	Immunology

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481