

## Anti-IL-1 beta/IL-1F2 Antibody (3V937)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3V937
Purification:	Protein G purified

### Applications

Application:	ELISA
--------------	-------

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human Interleukin-1 beta
Antigen Species:	Human
Uniprot ID:	P01584
Synonyms:	Interleukin-1 beta;IL1F2;IL-1 beta;IL1B;Interleukin-1 $\beta$ ;IL-1B;IL-1 $\beta$ /IL-1F2;IL-1BETA;IL-1 $\beta$ ;IL-1 $\beta$
Biology Area:	Immunology

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481