

## Anti-KAAG1 Antibody (7Y784)

### Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7Y784
Purification:	Affinity-chromatography

### Applications

Verified Activity:	The Binding Activity of Human KAAG1 with Anti-KAAG1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human KAAG1 at 2 µg/mL can bind Anti-KAAG1 recombinant antibody. The EC50 is 2.040-2.284 ng/mL.
Application:	ELISA

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human KAAG1 Protein
Antigen Species:	Human
Uniprot ID:	Q9UBP8
Biology Area:	Immunology

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481