

## Anti-Kanamycin Antibody (3J6)

### Product Details

Ig Type:	IgG1
Conjugation:	Unconjugated
Clone:	3J6
Purification:	Protein G purified

### Applications

Application:	ELISA
--------------	-------

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Kanamycin-BSA conjugated
Synonyms:	Amkin;enterokan;biograncetin
Biology Area:	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481