

Anti-KRT19 Antibody (7H611)

Product Details

Ig Type:	IgG2a
Reactivity:	Human, Mouse
Conjugation:	Unconjugated
Clone:	7H611
Purification:	Protein G purified

Applications

Verified Activity:	<ol style="list-style-type: none">Western Blot<ul style="list-style-type: none">-All lanes: Mouse anti-Human Cyfra21-1(KS19.1) monoclonal antibody at 4ug/ml-Lane 1: MCF-7 cell lysate-Lane 2: MDA-MB-231 cell lysate-Lane 3: HepG2 cell lysate-Secondary HRP labeled Goat polyclonal to Mouse IgG at 1/3000 dilution-Predicted band size : 30 kDa-Observed band size : 32 kDa-Additional bands at: 50 kDa <ol style="list-style-type: none">Immunohistochemical of paraffin-embedded human prostate tissue using TMAH-00675 at dilution of 1:200Immunohistochemical of paraffin-embedded human esophageal tissue using TMAH-00675 at dilution of 1:200
Application:	ELISA,IHC,WB

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Cyfra21-1 (KS19.1)
Antigen Species:	Human
Uniprot ID:	P08727
Biology Area:	Signal Transduction

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481