

## Anti-Lipocalin-2/LCN2 Antibody (5Z214)

### Product Details

Ig Type:	IgG2b
Reactivity:	Human, Mouse
Conjugation:	Unconjugated
Clone:	5Z214
Purification:	Protein G purified

### Applications

Verified Activity:	1. Western Blot
	-All lanes: Mouse anti-human Neutrophil gelatinase-associated lipocalin monoclonal antibody at 1µg/ml
	-Lane 1: NGAL transfected 293 cell lysate
	-Predicted band size : 22 kDa
	-Observed band size : 25 kDa
	2. Immunohistochemical of paraffin-embedded human spleen organization using TMAH-00685 at dilution of 1:200
	3. Immunohistochemical of paraffin-embedded human prostate tissue using TMAH-00685 at dilution of 1:200
Application:	ELISA, WB

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human Neutrophil gelatinase-associated lipocalin Protein
Antigen Species:	Human
Uniprot ID:	P80188
Synonyms:	Lipocalin-2;MSFI;p25;Oncogene 24p3;Siderocalin LCN2;NGAL
Biology Area:	Immunology

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481