

Anti-LGALS3 Antibody (5C881)

Product Details

Ig Type:	IgG2b
Reactivity:	Human, Mouse
Conjugation:	Unconjugated
Clone:	5C881
Purification:	Protein G purified

Applications

1. Immunohistochemical of paraffin-embedded human ovarian cancer using TMAH-00692 at dilution of 1:100
2. Immunohistochemical of paraffin-embedded human small intestinal tissue using TMAH-00692 at dilution of 1:100
3. Western Blot

Verified Activity:	-Positive WB detected in: MCF-7 whole cell lysate, Hela whole cell lysate, NIH/3T3 whole cell lysate -All lanes: LGALS3 antibody at 1:2000 -Secondary: Goat polyclonal to Mouse IgG at 1/50000 dilution -Predicted band size: 27 kDa -Observed band size: 27 kDa
--------------------	--

Application: ELISA,IHC,WB

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Galectin 3(GAL-3)
Antigen Species:	Human
Uniprot ID:	P17931
Synonyms:	Gal-3;CBP35;GALIG;MAC2GAL3;Galectin3;Lectin L-29;Mac-2;L29;LGALS3;L31;GAL3;AGE-R3;galactin-3;LGALS2;GALBPCBP35
Biology Area:	Neuroscience

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481