

Anti-LMO2 Antibody (3P24)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3P24
Purification:	Affinity-chromatography

Applications

	Western Blot
	-Positive WB detected in: U937 whole cell lysate, U87 whole cell lysate, PC3 whole cell lysate, MCF-7 whole cell lysate
Verified Activity:	-All lanes: LMO2 antibody at 1:2000
	-Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution
	-Predicted band size: 19, 26, 6 kDa
	-Observed band size: 20-25, 90-100 kDa
Application:	ELISA,WB
Recommended	WB:1:500-1:2000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human LMO2
Antigen Species:	Human
Gene ID:	4005
Uniprot ID:	P25791
Synonyms:	RBTN2;TTG2;RHOM2;RBTNL1;T-cell translocation protein 2;Rhombotin-2;Cysteine-rich protein TTG-2;LIM domain only protein 2;LMO-2
Biology Area:	Cardiovascular

Research Background

Acts with TAL1/SCL to regulate red blood cell development. Also acts with LDB1 to maintain erythroid precursors in an immature state.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481