

Anti-MMBP Antibody (2K358)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Escherichia coli
Conjugation:	Unconjugated
Clone:	2K358
Purification:	Affinity-chromatography

Applications

	Western Blot
	-Positive WB detected in Escherichia coli lysate
Verified Activity:	-All lanes: malE antibody at 1: 1000
	-Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution
	-Predicted band size: 44 kDa
	-Observed band size: 44 kDa
Application:	ELISA,WB

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Escherichia coli malE
Antigen Species:	Escherichia coli (strain K12)
Uniprot ID:	P0AEX9
Biology Area:	Tags & Cell Markers

Research Background

Part of the ABC transporter complex MalEFGK involved in maltose/maltodextrin import. Binds maltose and higher maltodextrins such as maltotriose.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481