

## Anti-MAP3K3 Antibody (5S171)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5S171
Purification:	Affinity-chromatography

### Applications

Application:	ELISA
--------------	-------

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthetic peptide: Human MAP3K3
Antigen Species:	Human
Gene ID:	4215
Uniprot ID:	Q99759
Synonyms:	Mitogen activated protein kinase kinase kinase 3;MEKK3;MEK kinase 3;MEKK 3;MAP/ERK kinase kinase 3;M3K3_HUMAN;MAPK/ERK kinase kinase 3;MAPKKK3
Biology Area:	Signal Transduction

### Research Background

Component of a protein kinase signal transduction cascade. Mediates activation of the NF-kappa-B, AP1 and DDIT3 transcriptional regulators.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481