

Anti-Mucin-16/MUC16 Antibody (9B20)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9B20
Purification:	Affinity-chromatography

Applications

Verified Activity:	IHC image of TMAH-00788 diluted at 1:100 and staining in paraffin-embedded human endometrial cancer performed on a Leica Bond TM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a Goat anti-rabbit IgG polymer labeled by HRP and visualized using 0.05% DAB.
Application:	ELISA, IHC
Recommended	IHC:1:50-1:200.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human CA-125 (MUC16)
Antigen Species:	Human
Gene ID:	94025
Uniprot ID:	Q8WXI7
Biology Area:	Cancer, Tags & Cell Markers, Signal transduction

Research Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481