

Anti-MYCBP Antibody (3Q306)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3Q306
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded Human testis tissue using TMAH-00794 at dilution of 1:200.
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human MYCBP Protein
Antigen Species:	Human
Gene ID:	26292
Uniprot ID:	Q99417
Synonyms:	AMY1;Associate of Myc 1;C myc binding protein;AMY 1;Cmyc binding protein;AMY-1
Biology Area:	Epigenetics and Nuclear Signaling

Research Background

May control the transcriptional activity of MYC. Stimulates the activation of E box-dependent transcription by MYC.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481