

Anti-Neomycin Antibody (5V885)

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Clone:	5V885
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide
Synonyms:	CCRIS 5462;HSDB 3242;Caswell No. 595;Myacyne.;Fradiomycinum;EPA Pesticide Chemical Code 006303;Fradiomycin
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481