

Anti-Ofloxacin Antibody (6Z515)

Product Details

Ig Type: IgG2b
Conjugation: Unconjugated
Clone: 6Z515
Purification: Protein G purified

Applications

Application: ELISA

Properties

Purity: >95%
Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Ofloxacin-BSA conjugated
Synonyms: orocin;Oflocee;oflocin;oflocet;tarivid
Biology Area: Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481