

Anti-Pepsinogen C/PGC Antibody (5U637)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5U637
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human gastric cancer tissue using TMAH-00878 at dilution of 1:200
Application:	ELISA,IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human PGC
Antigen Species:	Human
Gene ID:	5225
Uniprot ID:	P20142
Synonyms:	PGII;progastricsin (pepsinogen C);PEPC
Biology Area:	Signal Transduction

Research Background

Hydrolyzes a variety of proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481