

Anti-PLEKHA1 Antibody (8U160)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8U160
Purification:	Affinity-chromatography

Applications

	Western Blot
	-Positive WB detected in: Hela whole cell lysate, A549 whole cell lysate, PC-3 whole cell lysate, MCF-7 whole cell lysate
Verified Activity:	-All lanes: PLEKHA1 antibody at 1:1000
	-Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution
	-Predicted band size: 46, 39 kDa
	-Observed band size: 36-55 kDa
Application:	ELISA,WB
Recommended	WB:1:500-1:2000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human PLEKHA1
Antigen Species:	Human
Gene ID:	59338
Uniprot ID:	Q9HB21
Synonyms:	PH domain-containing family A member 1;Pleckstrin homology domain-containing family A member 1;TAPP-1;Tandem PH domain-containing protein 1;TAPP1
Biology Area:	Signal transduction

Research Background

Binds specifically to phosphatidylinositol 3,4-diphosphate (PtdIns3,4P2), but not to other phosphoinositides. May recruit other proteins to the plasma membrane.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481