

Anti-PRSS2 Antibody (3Y816)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3Y816
Purification:	Affinity-chromatography

Applications

Verified Activity:	Immunofluorescence staining of Jurkat cell with TMAH-01010 at 1:50, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 562-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).
Application:	ELISA, IF
Recommended	IF:1:50-1:200.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human PRSS2
Antigen Species:	Human
Gene ID:	5645
Uniprot ID:	P07478
Biology Area:	Metabolism, Signal transduction

Research Background

In the ileum, may be involved in defensin processing, including DEFA5.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481