

Anti-PXN Antibody (2F835)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2F835
Purification:	Affinity-chromatography

Applications

Application:	ELISA
--------------	-------

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human PXN
Antigen Species:	Human
Gene ID:	5829
Uniprot ID:	P49023
Synonyms:	Paxillin;PAX1;Paired box protein Pax 1;FLJ23042;FLJ16691;PAXI_HUMAN;PAX 1;Paxillin alpha
Biology Area:	Signal Transduction

Research Background

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481