

Anti-RRM1 Antibody (4N328)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4N328
Purification:	Protein G purified

Applications

Verified Activity:	1. Immunohistochemistry of paraffin-embedded human lung cancer using TMAH-01056 in 30µg/ml dilute concentrations.
	2. Immunohistochemistry of paraffin-embedded human breast cancer using TMAH-01056 in 30µg/ml dilute concentrations.
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human RRM1 Protein
Antigen Species:	Human
Gene ID:	6240
Uniprot ID:	P23921
Synonyms:	ribonucleotide reductase M1;RR1;RIR1;R1
Biology Area:	Epigenetics and Nuclear Signaling

Research Background

Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the corresponding ribonucleotides.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481