

Anti-SDC1 Antibody (8E385)

Product Details

Ig Type:	IgG1
Reactivity:	Rat
Conjugation:	Unconjugated
Clone:	8E385
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Rat Syndecan-1 Protein
Antigen Species:	Rat
Gene ID:	25216
Uniprot ID:	P26260
Synonyms:	syndecan 1
Biology Area:	Neuroscience

Research Background

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix. Regulates exosome biogenesis in concert with SDCBP and PDCD6IP.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481