

Anti-SDC1 Antibody (8J461)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8J461
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SDC1 Protein
Antigen Species:	Human
Gene ID:	6382
Uniprot ID:	P18827
Synonyms:	SDC;syndecan 1;syndecan;SYND1;Syndecan-1;CD138
Biology Area:	Neuroscience

Research Background

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix. Regulates exosome biogenesis in concert with SDCBP and PDCD6IP.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481