

Anti-SLC39A6 Antibody (8F657)

Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8F657
Purification:	Affinity-chromatography

Applications

Verified Activity:	<p>1. The Binding Activity of SLC39A6 with anti-SLC39A6 antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SLC39A6 at 1 µg/ml can bind Anti-SLC39A6-1 recombinant antibody, the EC50 is 0.6873-0.9010 ng/mL.</p> <p>2. The Binding Activity of Cynomolgus SLC39A6 with Anti-SLC39A6 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Cynomolgus SLC39A6 at 2 µg/mL can bind Anti-SLC39A6 recombinant antibody. The EC50 is 0.8290-0.9864</p>
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SLC39A6 Protein
Antigen Species:	Human
Gene ID:	25800
Uniprot ID:	Q13433
Synonyms:	ZIP6;LIV-1;LIV1
Biology Area:	Immunology

Research Background

May act as a zinc-influx transporter.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481