

Anti-SN38 Antibody (7U210)

Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7U210
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Labetuzumab govitecan (Synonyms: IMMU 130; hMN-14-SN-38) with Anti-SN38 antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Labetuzumab govitecan (Synonyms: IMMU 130; hMN-14-SN-38) at 2 µg/mL can bind Anti-SN38 antibody, the EC50 is 3.326 to 4.721 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	SN38-BSA
Synonyms:	7-Ethyl-10-hydroxycamptothecin
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481