

Anti-Osteopontin Antibody (6U756)

Product Details

| | |
|---------------|--------------------|
| Ig Type: | IgG2a |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 6U756 |
| Purification: | Protein G purified |

Applications

1. Immunohistochemical of paraffin-embedded human colon cancer using TMAH-01128 at dilution of 1:200
2. Immunohistochemical of paraffin-embedded human Kidney tissue using TMAH-01128 at dilution of 1:200
3. Western Blot

Verified Activity:

- All lanes Mouse anti human OPN monoclonal antibody at at 2µg/ml
- Lane 1:mouse kidney tissue lysate
- Secondary Goat polyclonal to Mouse IgG at 1/5000 dilution
- Predicted band size : 33,36,34kDa
- Observed band size: 50kDa

Application: ELISA, WB, IHC

Properties

Purity: >95%

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

| | |
|------------------|---|
| Immunogen: | Recombinant Protein: Human Osteopontin |
| Antigen Species: | Human |
| Uniprot ID: | P10451 |
| Synonyms: | ETA-1;OPN;BSPI;BNSP;secreted phosphoprotein 1 |
| Biology Area: | Cardiovascular |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481