

Anti-STAT3 Antibody (1E44)

Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1E44
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Human STAT3 with Anti-STAT3 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human STAT3 at 2 µg/mL can bind Anti-STAT3 recombinant antibody, the EC50 is 41.31-75.59 ng/mL.
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human STAT3 Protein
Antigen Species:	Human
Gene ID:	6774
Uniprot ID:	P40763
Synonyms:	signal transducer and activator of transcription 3 (acute-phase response factor);APRF;ADMIO; HIES
Biology Area:	Epigenetics and Nuclear Signaling, Cancer, Cardiovascular, Developmental biology, Signal transduction, Stem cells

Research Background

Signal transducer and transcription activator that mediates cellular responses to interleukins, KITLG/SCF, LEP and other growth factors. Once activated, recruits coactivators, such as NCOA1 or MED1, to the promoter region of the target gene. May mediate cellular responses to activated FGFR1, FGFR2, FGFR3 and FGFR4. Upon activation of IL6ST/gp130 signaling by interleukin-6 (IL6), binds to the IL6-responsive elements identified in the promoters of various acute-phase protein genes. Activated by IL31 through IL31RA. Acts as a regulator of inflammatory response by regulating differentiation of naive CD4(+) T-cells into T-helper Th17 or regulatory T-cells (Treg): deacetylation and oxidation of lysine residues by LOXL3, leads to disrupt STAT3 dimerization and inhibit its transcription activity. Involved in cell cycle regulation by inducing the expression of key genes for the progression from G1 to S phase, such as CCND1. Mediates the effects of LEP on melanocortin production, body energy homeostasis and lactation. May play an apoptotic role by transactivating BIRC5 expression under LEP activation. Cytoplasmic STAT3 represses macroautophagy by inhibiting EIF2AK2/PKR activity. Plays a crucial role in basal beta cell functions, such as regulation of insulin secretion.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481