

Anti-STK38L Antibody (7Z651)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7Z651
Purification:	Protein G purified

Applications

Western Blot

- Positive WB detected in: SKT38L antibody at 1:1000
- Lane 1: HepG2 whole cell lysate
- Lane 2: 293T whole cell lysate tissue
- Lane 3: Jurkat whole cell lysate
- Lane 4: Raw264.7 whole cell lysate
- Lane 5: K562 whole cell lysate
- Secondary: Goat polyclonal to Mouse IgG at 1/20000 dilution
- Predicted band size: 55, 44 KDa
- Observed band size: 55 KDa
- Exposure time: 5min

Verified Activity:

Application: ELISA, WB

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Serine/threonine-Protein kinase 38-like Protein (212-464AA)
Antigen Species:	Human
Gene ID:	23012
Uniprot ID:	Q9Y2H1
Biology Area:	Others

Research Background

Involved in the regulation of structural processes in differentiating and mature neuronal cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481