

Anti-Sulfadimidine Antibody (1M267)

Product Details

Ig Type:	IgG2a
Conjugation:	Unconjugated
Clone:	1M267
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Sulfadimidine-BSA conjugated
Synonyms:	4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulf-onamide;Sulfamethazine,4,6-Dimethylsulfadiazine;S-Mez;SM2
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481