

Anti-Sulfametoxydiazine Antibody (1Q847)

Product Details

Ig Type:	IgG1
Conjugation:	Unconjugated
Clone:	1Q847
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Sulfametoxydiazine-BSA conjugated
Synonyms:	SMD;durenat;bayrena;aristebon;5-Methoxysulfadiazine;sulfameter
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481