

## Anti-Tissue Factor Antibody (5C478)

### Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5C478
Purification:	Affinity-chromatography

### Applications

Verified Activity:	<p>1. The Binding Activity of Human TF with Anti-TF Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TF at 2 µg/mL can bind Anti-TF recombinant antibody. The EC50 is 1.434-1.635 ng/mL.</p> <p>2. The Binding Activity of Cynomolgus TF with Anti-TF Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Cynomolgus TF at 2 µg/mL can bind Anti-TF recombinant antibody. The EC50 is 0.4495-0.5538 ng/mL.</p>
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human TF Protein
Antigen Species:	Human
Uniprot ID:	P13726
Biology Area:	Cardiovascular, Metabolism

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481