

## Anti-TFF3 Antibody (7Z596)

### Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7Z596
Purification:	Protein G purified

### Applications

Verified Activity:	1. Immunohistochemical of paraffin-embedded human colon cancer using TMAH-01163 at dilution of 1:50
	2. Immunofluorescent analysis of MCF-7 cells using TMAH-01163 at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)
Application:	ELISA,IF,IHC

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human Trefoil factor 3
Antigen Species:	Human
Gene ID:	7033
Uniprot ID:	Q07654
Biology Area:	Signal Transduction

### Research Background

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481