

Anti-TM4SF1 Antibody (2G843)

Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2G843
Purification:	Affinity-chromatography

Applications

Verified Activity:	<p>1. The Binding Activity of Human TM4SF1 with Anti-TM4SF1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TM4SF1 at 5 µg/mL can bind Anti-TM4SF1 recombinant antibody. the EC50 is 4.079-4.472 ng/mL.</p> <p>2. Untransfected CT26 cells (green line) and transfected Human TM4SF1*10⁶ cells, washed and then followed by APC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.</p>
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human TM4SF1 Protein
Antigen Species:	Human
Gene ID:	4071
Uniprot ID:	P30408
Synonyms:	H-L6;TAAL6;M3S1;L6
Biology Area:	Cancer

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481