

Anti-TNNI3 Antibody (7G122)

Product Details

Ig Type:	IgG2b
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7G122
Purification:	Protein G purified

Applications

Application:	ELISA
--------------	-------

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Cardiac Troponin I Protein
Antigen Species:	Human
Gene ID:	7137
Uniprot ID:	P19429
Biology Area:	Signal Transduction

Research Background

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481