

Anti-TSLP Antibody (40865)

Product Details

Ig Type:	Human IgG2
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	40865
Purification:	Affinity-chromatography

Applications

Application:	ELISA
--------------	-------

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human TSLP Protein
Antigen Species:	Human
Gene ID:	85480
Uniprot ID:	Q969D9
Synonyms:	thymic stromal lymphopoietin;TSLP
Biology Area:	Inflammation, Cancer

Research Background

Cytokine that induces the release of T-cell-attracting chemokines from monocytes and, in particular, enhances the maturation of CD11c(+) dendritic cells. Can induce allergic inflammation by directly activating mast cells. May act as an antimicrobial peptide in the oral cavity and on the skin.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481