

## Anti-TSLP Antibody (6F714)

### Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6F714
Purification:	Affinity-chromatography

### Applications

Verified Activity:	The Binding Activity of Human TSLP(R127A,R130A) with Anti-TSLP recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized human TSLP(R127A, R130A) at 2 µg/mL can bind Anti-TSLP recombinant antibody. The EC50 is 52.55-58.67 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human TSLP Protein
Antigen Species:	Human
Uniprot ID:	Q969D9
Synonyms:	thymic stromal lymphopoietin;TSLP
Biology Area:	Immunology

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481