

Anti-TSPAN8 Antibody (6V251)

Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6V251
Purification:	Affinity-chromatography

Applications

Verified Activity:	<p>1. The Binding Activity of TSPAN8 with Anti-TSPAN8 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TSPAN8 at 5 µg/mL can bind Anti-TSPAN8 recombinant antibody. The EC50 is 2.261-2.623 ng/mL. The VLPs is negative control.</p> <p>2. Untransfected MC38 cells (green line) and transfected Human TSPAN8 MC38 stable cells (red line) were stained with anti-TSPAN8 antibody (TMAH-01221) (2µg/1*10⁶ cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.</p> <p>3. The purity of TSPAN8 was greater than 95% as determined by SEC-HPLC</p>
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human TSPAN8 Protein
Antigen Species:	Human
Gene ID:	7103
Uniprot ID:	P19075
Biology Area:	Cancer, Tags & Cell Markers, Signal transduction

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481