

Anti-ULBP-1 Antibody (2Q93)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2Q93
Purification:	Affinity-chromatography

Applications

	Western Blot
	-Positive WB detected in: Hela whole cell lysate, MCF-7 whole cell lysate, Raji whole cell lysate, HEK293 whole cell lysate
Verified Activity:	-All lanes: ULBP1 antibody at 1:1000
	-Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution
	-Predicted band size: 28 kDa
	-Observed band size: 36-50 kDa
Application:	ELISA, WB
Recommended	WB:1:500-1:2000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human ULBP1
Antigen Species:	Human
Gene ID:	80329
Uniprot ID:	Q9BZM6
Synonyms:	RAET11;UL16 binding protein 1
Biology Area:	Immunology

Research Background

Binds and activates the KLRK1/NKG2D receptor, mediating natural killer cell cytotoxicity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481