

Anti-VSIG4 Antibody (3W528)

Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3W528
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of VSIG4 with anti-VSIG4 antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human VSIG4 at 2 µg/ml can bind Anti-VSIG4 recombinant antibody, the EC50 is 51.14-68.73 ng/mL.
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human VSIG4 Protein
Antigen Species:	Human
Gene ID:	11326
Uniprot ID:	Q9Y279
Synonyms:	CRlg;V-set and immunoglobulin domain containing 4;Z39IG
Biology Area:	Immunology

Research Background

Phagocytic receptor, strong negative regulator of T-cell proliferation and IL2 production. Potent inhibitor of the alternative complement pathway convertases.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481